



## **MERITS INTERCEPTOR SCOOTER**

**\$8,194.00**



**SKU: N/A**

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

This large sporty streamline design scooter has full suspension with shock absorbers for a more comfortable ride. Your battery gauge, temperature, accumulated usage hours and maintenance messages are all displayed on the LCD digital dashboard. This scooter also has a handy storage space at the tiller and a comfortable seating system.

### **FEATURES**

- Brand new streamline design with a sporty feel
- Fully adjustable suspension system that delivers a smooth and comfortable ride
- Modern LED lighting system delivers brilliant light with less power consumption
- Digital LCD dashboard display including a time display, battery gauge, temperature, maintenance messages, speed display and odometer
- Extra storage space beneath tiller
- Built-in half speed switch which operates automatically when turning to provide extra safety
- Merits 4PT (4 pole technology) and large controller, delivering more power for climbing steep inclines
- New easy to access 12 volt socket

### **SPECIFICATIONS**

- Overall Length: 154 cm
- Overall Width: 72 cm
- Overall Height: 145.5 cm
- Seat height from ground: 60-68 cm

- Seat height from deck: 43-50 cm
- Wheel base: 113 cm
- Total Capacity: 204 kg (450 lbs)
- Max Speed: 10-15 km
- Max Range per charge: 40-60 km
- Turn Radius: 180 cm
- Ground clearance: 13 cm
- Motor: DC 24 V, 650W (rating) / 220 W (max)
- Controller: PG S-drive 200A
- Batteries: 12V , 60-100AH x 2
- Charger: 8A off-board
- Gradient: 10 degrees
- Tires: Pneumatic
  - front 14"
  - rear 16"
  - Brakes: Intelligent, regenerative, electromagnetic brakes
- Total Weight with Batteries: 148 kg (325 lbs)

Buy This Merits Interceptor Scooter online at MediMart using Afterpay or [NDIS method](#) today! Expert [mobility equipment suppliers](#) in Australia.