



## MERITS MAVERICK 12 REHAB POWER CHAIR



**SKU:** P327A-3-39905210

## PRODUCT DESCRIPTION

### Features

- Transit Occupied Certified
- Mid-wheel drive design offers outstanding manoeuvrability in tight, compact spaces and the superior adjustable suspension design delivers a smooth and stable ride both indoors and outdoors, even over tough terrain
- Patented obstacle climbing suspension that works
- Dual in-line motors for enhanced efficiency, torque, range, and performance
- Height & depth adjustable armrests
- Height adjustable four post seats, with width & depth adjustment (R601)
- After market cushions & backrests can be easily fitted to the seat
- 204 kg weight capacity
- Available in a number of vibrant new colours – Red, Silver, Pearl Gold, Charcoal, NV Green, Kite Blue
- Crash tested tie down points
- Other options available are head rest, rear attendant controls, optional leg rests

### Specifications

- Overall Dimension : A x B x C 90x 69 x 100cm / 35 x 27 x 39in
- Seat Dimension : D x E x F 46-56 x 46-56 x 46cm / 18-22 x 18-22 x 18in
- Seat Height from Ground : G 44-49cm / 17-19in

- Seat Height from Deck : H 32-37cm / 12.5-14.5in
- Max Speed 8kph / 5mph
- Range 32km / 20mi.
- Turn Radius 53cm / 21"
- Ground clearance 5.6cm / 2.2"
- Motor DC24V, 200W(Rating) / 575W(Max)
- Controller PG VR2 50A
- Battery 12V / 55-62A/H x 2pcs
- Charger 5A off-board Gradient 8°
- Caster front: 6" x 2" PU tire ; rear: 6" x 2" solid tire
- Drive Wheel 12" (4.00-5) pneumatic tire
- Brake intelligent, regenerative, electromagnetic brakes