




MERITS P312A DUALER POWER CHAIR



SKU: P31231WRMA

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Easily turn around in small places with the Dualer Power Chair with an elevating seat. The turning radius is 21" and easy to control by using the swing-away joystick. This versatile power chair reaches the top speed of 8 kmph, has a weight capacity of 136 kilos, and has remarkable handling.

FEATURES

- Easy to convert from RWD to FWD, 2-in-1 power chair, ideal chair for indoor as well as outdoor use.
- Unique compact size for maximum manoeuvrability in narrow space and walkway.
- Height adjustable flip-up armrests can provide seat width adjusted.
- Start with RWD and move to FWD by rotating the chair 180°, to enjoy direct-drive, direct-response feeling.
- Turning radius is only 53cm/ 21".
- Low seat height of 18" provides excellent stability.
- Power seat elevator available.

SPECIFICATIONS

Driving Range	28 km / 17.5 mil
Max Speed	8 km/h , 5 mph

Battery Type	12V / U1 x 2pcs
Motor	DC 24V, 1240W (Rating) / 950W (Max)
Max Climbing Angle	10°
Turning Radius	53 cm / 21"
Ground Clearance	5 cm / 2"
Controller	PG VS2 60
Charger	5A off-board
Weight Capacity	136 kg / 300 lb
Total Weight	85 kg / 187 lb
Front wheel	8" (200 x 50) PU tire
Drive Wheel	10" (260 x 85) foam filled tire
Rear wheel	10" (260 x 85) Foam filled tire

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au