




MERITS WRANGLER SCOOTER



SKU: S840DXG-BB-H

PRODUCT DESCRIPTION

Purpose & Benefits of Merits Wrangler Scooter

The **Merits Wrangler Scooter** is a robust, high-performance mobility scooter built for users who require maximum power, comfort, and durability for demanding outdoor conditions. Designed with heavy-duty suspension, large pneumatic tyres, and a high ground clearance, the Wrangler is ideal for tackling rough terrain, long-distance travel, and uneven outdoor paths with ease.

With a weight capacity of 181 kg, this scooter provides enhanced stability and user confidence. The ergonomic seat, full lighting package, and advanced safety features ensure a secure and comfortable ride, making it a reliable choice for active users who need a dependable mobility solution for everyday outdoor use.

Key Features

- Heavy-duty performance for rugged outdoor terrain
- High ground clearance and large pneumatic tyres for a smoother ride
- Full suspension system for improved comfort and shock absorption
- Comfortable, adjustable seat with headrest and armrests
- Ergonomic tiller with easy-to-use controls

- Front and rear lighting system with indicators and hazard lights
- 181 kg weight capacity for added strength and support
- Ideal for users who need long-range mobility and off-road capability

Specifications

| | |
|------------------------------|---------------------|
| Overall Length | 1290 mm |
| Overall Width | 610 mm |
| Overall Height | 1280 mm |
| Wheelbase Length | 950 mm |
| Seat Width | 460 mm |
| Seat Depth | 400 mm |
| Min Floor To Seat Height | 430 mm |
| Max Floor To Seat Height | 490 mm |
| Min Seat Height from Deck | 640 mm |
| Max Seat Height from Deck | 710 mm |
| Turning Radius | 1500 mm |
| Backrest Height | 510 mm |
| minrangeKM | 30 km |
| maxrangeKM | 40 km |
| Max Speed | 10 km/h |
| Front Wheel | 12" pneumatic |
| Rear Wheel | 12" pneumatic |
| Gradient | 8 ° |
| Ground Clearance | 120 mm |
| Battery | 12V/50A/H x 2pcs Ah |
| Motor Size | 400W/1750W hp |
| Controller | Rhino 2 120A |
| Charger | 4Amp A |
| Product Weight | 100 kg |
| Total Weight Without Battery | 72 kg |
| Max User Weight | 159 kg |
| Brand | Merits |

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au