




MERITS WRANGLER SCOOTER

~~\$5,104.00~~ Original price was:

~~\$5,104.00.~~ \$4,593.24 Current price is: \$4,593.24.



SKU: S840DXG-BB-H

PRODUCT DESCRIPTION

Purpose & Benefits of Merits Wrangler Scooter

The **Merits Wrangler Scooter** is a robust, high-performance mobility scooter built for users who require maximum power, comfort, and durability for demanding outdoor conditions. Designed with heavy-duty suspension, large pneumatic tyres, and a high ground clearance, the Wrangler is ideal for tackling rough terrain, long-distance travel, and uneven outdoor paths with ease.

With a weight capacity of 181 kg, this scooter provides enhanced stability and user confidence. The ergonomic seat, full lighting package, and advanced safety features ensure a secure and comfortable ride, making it a reliable choice for active users who need a dependable mobility solution for everyday outdoor use.

Key Features

- Heavy-duty performance for rugged outdoor terrain
- High ground clearance and large pneumatic tyres for a smoother ride
- Full suspension system for improved comfort and shock absorption
- Comfortable, adjustable seat with headrest and armrests
- Ergonomic tiller with easy-to-use controls

- Front and rear lighting system with indicators and hazard lights
- 181 kg weight capacity for added strength and support
- Ideal for users who need long-range mobility and off-road capability

Specifications

Overall Length	1290 mm
Overall Width	610 mm
Overall Height	1280 mm
Wheelbase Length	950 mm
Seat Width	460 mm
Seat Depth	400 mm
Min Floor To Seat Height	430 mm
Max Floor To Seat Height	490 mm
Min Seat Height from Deck	640 mm
Max Seat Height from Deck	710 mm
Turning Radius	1500 mm
Backrest Height	510 mm
minrangeKM	30 km
maxrangeKM	40 km
Max Speed	10 km/h
Front Wheel	12" pneumatic
Rear Wheel	12" pneumatic
Gradient	8 °
Ground Clearance	120 mm
Battery	12V/50A/H x 2pcs Ah
Motor Size	400W/1750W hp
Controller	Rhino 2 120A
Charger	4Amp A
Product Weight	100 kg
Total Weight Without Battery	72 kg
Max User Weight	159 kg
Brand	Merits