




### **MERITS YOGA FOLDABLE SCOOTER**

**~~\$4,401.00~~ Original price was:**

**~~\$4,401.00.~~ \$3,961.28 Current price is: \$3,961.28.**



**SKU: S542**

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits of Merits Yoga Foldable Scooter**

The **Merits Yoga Foldable Scooter** is a lightweight and innovative mobility solution designed for users who need portability without compromising comfort or performance. Its unique folding frame and detachable battery system make it perfect for travel, storage, and everyday convenience — whether at home, in a car boot, or on holiday.

With a stylish magnesium alloy frame, the Yoga is built to be both durable and ultra-lightweight. It folds quickly and effortlessly, offering independence to users who want a compact mobility option for urban environments, shopping centres, or transport hubs. Despite its portability, the Yoga still provides a smooth, stable ride with rear suspension and a comfortable seating design.

### **Key Features:**

- Unique folding frame for easy transport and storage
- Ultra-lightweight magnesium alloy construction
- Detachable lithium battery for convenient charging
- Rear suspension for improved ride comfort
- Compact design ideal for travel and tight spaces

- Comfortable padded seat and easy-to-use controls
- Suitable for indoor and light outdoor use
- Airline-friendly battery (check with your carrier)

## **Specifications:**

- Overall Dimension – Unfolded : 98 × 43.5 × 102.5 cm
- Overall Dimension – Folded : 52 × 43.5 × 64 cm
- Seat Dimension : 40.5 × 35 × 31.5 cm
- Seat Height from Ground : 59.5 cm
- Seat Height from Deck : 39 cm
- Wheelbase : 78.5 cm
- Max Speed 6kph
- Range 18 km
- Turn Radius 96 cm
- Ground clearance 9.5 cm
- Motor DC 24V, 80W (Rating)/ 340W (Max)
- Controller PG S-Drive 45A
- Battery Li-Ion Battery 24V/11.5A/H
- Charger 2A off-board Gradient 6°
- Front Wheel 6”(153×40) foam filled tire
- Rear Wheel 7”(178×45) foam filled tire
- Brake Electromagnetic brake