



MERITS YOGA TRAVEL SCOOTER



SKU: N/A

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Merits Yoga is a lightweight electric travel scooter that can be folded up in a matter of seconds and wheeled along the ground or put into a suitcase and transported.

Great for active users who like to travel for use on relatively flat indoor and outdoor surfaces, it will get you around the shops with ease. It even features automatic reduced speed on turning!

Maximum speed of 6 km/h and a range per charge of up to 18 km, the Merits Yoga provides a comfortable and convenient ride.

This space-saving scooter folds down quickly and easily and comes for a limited time with a **FREE suitcase** to transport it in.

FEATURES

- Strong, light-weight alloy aluminum frame designed for easy folding and carrying
- Small, foldable scooter for easy storage and transportation
- Lightweight 24v Li-ion battery with quick charge 2a off board charger
- Height and angle adjustable steering tiller to fit individual users
- Ideal for both indoor and outdoor use
- Patented Twin-front-wheels w/spring suspension system offers excellent stability
- Armrest as standard.
- Comfortable padded seat and backrest
- Optional battery charging station, for charging the battery away from the scooter.

SPECIFICATIONS

Total Size A x B x C	980mm x 435mm x 1025mm
Total Folded Size A x B x C	520mm x 435mm x 640mm
Seat Size D x E x F	405mm x 350mm x 315mm
Seat Height from Ground G	595mm
Seat Height from Deck H	390mm
Wheel Base I	785mm
Maximum Speed	6kph
Range	18km
Weight with Battery, without Armrest	25.7kg
Battery Weight	1.7kg
SWL	113kg

Buy This Merits Yoga Travel Scooter online at MediMart using Afterpay or [NDIS method](#) today! Expert [mobility equipment suppliers](#) in Australia.

For further information go to [medimart.com.au](https://www.medimart.com.au) or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au