



## MERITS YOGA TRAVEL SCOOTER



SKU: N/A

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Merits Yoga is a lightweight electric travel scooter that can be folded up in a matter of seconds and wheeled along the ground or put into a suitcase and transported.

Great for active users who like to travel for use on relatively flat indoor and outdoor surfaces, it will get you around the shops with ease. It even features automatic reduced speed on turning!

Maximum speed of 6 km/h and a range per charge of up to 18 km, the Merits Yoga provides a comfortable and convenient ride.

This space-saving scooter folds down quickly and easily and comes for a limited time with a **FREE suitcase** to transport it in.

### FEATURES

- Strong, light-weight alloy aluminum frame designed for easy folding and carrying
- Small, foldable scooter for easy storage and transportation
- Lightweight 24v Li-ion battery with quick charge 2a off board charger
- Height and angle adjustable steering tiller to fit individual users
- Ideal for both indoor and outdoor use
- Patented Twin-front-wheels w/spring suspension system offers excellent stability
- Armrest as standard.
- Comfortable padded seat and backrest
- Optional battery charging station, for charging the battery away from the scooter.

## SPECIFICATIONS

<b>Total Size   A x B x C</b>	980mm x 435mm x 1025mm
<b>Total Folded Size   A x B x C</b>	520mm x 435mm x 640mm
<b>Seat Size   D x E x F</b>	405mm x 350mm x 315mm
<b>Seat Height from Ground   G</b>	595mm
<b>Seat Height from Deck   H</b>	390mm
<b>Wheel Base   I</b>	785mm
<b>Maximum Speed</b>	6kph
<b>Range</b>	18km
<b>Weight with Battery, without Armrest</b>	25.7kg
<b>Battery Weight</b>	1.7kg
<b>SWL</b>	113kg

Buy This Merits Yoga Travel Scooter online at MediMart using Afterpay or [NDIS method](#) today! Expert [mobility equipment suppliers](#) in Australia.