




MERTIS BARIATRIC FOLDING WALKER WITH WHEEL KIT



SKU: W144

PRODUCT DESCRIPTION

Purpose & Benefits

The **Merits W144 Walking Frame** is designed for bariatric users requiring stable and supportive mobility assistance. Its oversized frame offers added strength and durability while remaining easy to fold for convenient transport and storage.

Features

- Oversized frame suitable for bariatric users
- Strong and durable construction for increased stability
- Easily folds for convenient transport and storage
- Supports users up to 200 kg

Specifications

- **Overall Width:** 585 mm
- **Overall Height:** 820 – 1000 mm (adjustable)
- **Overall Length:** 550 mm
- **Unit Weight:** 3.5 kg

- **Load Capacity:** 200 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au