



### **MOLIFT AIR 205 MOTOR CEILING HOIST**



SKU: MLA205M

### PRODUCT DESCRIPTION

### **PURPOSE AND BENEFITS**

The Molift Air 205 is a new generation ceiling hoist/lift that lifts and transfers patients in a safe, smooth and comfortable way. This ceiling hoist is innovative, functional and well designed. The Molift Air 205 ceiling hoist is suitable for lifting and transferring from a seated position as well as gait retraining and standing.

It is available with a comprehensive range of accessories. The Molift Air has a quick release coupling that does not require special tools for installation. It also is light in weight, making it easy to install and remove. This ceiling hoist is available with a range of options such as remote control, In rail charging and weight scales.

It is suitable for ceiling hoist installations in home care, aged care facilities and hospitals.

Warranty on the Molift Air 205 is 2 years, as with all Molift Lifters.

Molift Batteries and chargers are under warranty for 6 months.

The Molift Air has a class leading lifting interval and is available in the following Safe Working Load options. Click on each of the sizes for more product information

# **Molift Air Weight Capacity Options**

- 205 kg
- 300 kg



### **FEATURES**

- Environmentally Friendly Ni MH batteries
- Automatic Service Notification Included
- In-Rail Charging
- Emergency Lowering and Emergency Stop
- Remote Control
- Propulsion/Traverse Motor
- Lifting Interval: 3,0 m/9 ft 10"
- Rated performance 50 hoists at 75 kg
- Hoisting Speed: 60mm/sec at 75kg load
- Emergency Lowering: Mechanical and Electrical
- Emergency-Stop: Reset Possible from Floor Level
- Lightweight Motor: 8.5 kg
- Modern Design
- High Weight Capacity
- Easy to Install
- Simple to Operate

For further information go to medimart.com.au or give us a call: 1300 633 063

## SA

# **Brighton Showroom**

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

#### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

## **Reynella Showroom**

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

revnella@medimart.com.au

#### VIC

# **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 319903 8763 3451

frankston@medimart.com.au

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168 03 7503 8354

monash@medimart.com.au



# NT

# **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 <a href="mailto:support@medimarnt.com.au">support@medimarnt.com.au</a>