

MOLIFT BARIATRIC 4-POINT SUSPENSION - PIG TAIL HOOK

\$1,950.00



SKU: 2510435



PRODUCT DESCRIPTION

Purpose & Benefits

The **Molift Bariatric 4-Point Suspension with Pig Tail Hook** is a durable and lightweight sling bar designed for safe and efficient bariatric patient transfers. Constructed from aluminium for low weight and easy handling, it offers a secure four-point lifting configuration to ensure optimal patient positioning and comfort during hoisting. Compatible with Molift Air 350 and Air 500 hoists, this suspension system enhances stability and ease of use for caregivers in bariatric care environments.

Features

- Lightweight aluminium construction for easy handling and reduced strain on caregivers
- Four-point suspension design ensures even weight distribution and patient comfort
- Pig tail hook system provides a secure and stable connection for slings
- Quick-release coupling allows fast and safe attachment or removal from hoists
- Compatible with Molift Air 350 and Air 500 ceiling hoists for bariatric lifting
- Designed for bariatric use, supporting up to 500 kg load capacity

Specifications

Specification	Details
Model	Molift Bariatric 4-Point Suspension – Pig Tail Hook
Brand	Molift
Compatibility	Molift Air 350 and Air 500 hoists
Material	Aluminium (lightweight, corrosion-resistant)
Load Capacity	Up to 500 kg
Overall Height	200 mm
Overall Width	550 mm
Overall Length	550 mm
Coupling Type	Quick release coupling for easy interchangeability