




MOLIFT DUO GANTRY FRAME



SKU: DuoGantry

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

A flexible gantry system

With a Molift freestanding unit, we offer a system that is easy to install without additional fixtures on ceiling or walls. It is ideally suited to temporary lift applications in hospitals, nursing homes, or private homes.

Lightweight construction

The Molift Rail System Duo has a lightweight aluminum construction and a functional design which makes the gantry system easy to install, dismantle and move elsewhere. In combination with Molift Air or Molift Nomad lift motor and a sling, Molift Rail System Duo solves all overhead lifting situations, whatever the environment. With a lifting capacity of 300 kg (660 lbs), it is ideal for lifting and transferring a person from or to; bed, floor, chair, wheelchair, toilet, etc.

Suitable for all bed and room sizes

This 2 post gantry system has a rail length between 200 cm (78,7") and 350 cm (137,8") and has 4 steps adjustable height up to 255 cm (100,4"), making this freestanding hoist unit suitable for all bed and room sizes.

SPECIFICATIONS

| | |
|-------------------------------|--------------------------|
| Available rail lengths | 200 - 250 - 300 - 350 cm |
| Height | 225-255 cm |

| | |
|-------------------------------|--|
| Depth | 120 cm |
| Maximum internal space | 350 x 240 cm |
| Safe Working Load SWL | 300 kg |
| Material | Anodized aluminum, plastic |
| Total Weight | 38-44,3 kg |
| Transport height | 170 cm |
| Weight | Post, complete: 14.8 kg Boom: 4.2 kg/m Foot: 3.8 kg Post, without foot: 11 kg |