



MOLIFT DUO GANTRY FRAME



SKU: DuoGantry

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

A flexible gantry system

With a Molift freestanding unit, we offer a system that is easy to install without additional fixtures on ceiling or walls. It is ideally suited to temporary lift applications in hospitals, nursing homes, or private homes.

Lightweight construction

The Molift Rail System Duo has a lightweight aluminum construction and a functional design which makes the gantry system easy to install, dismantle and move elsewhere. In combination with Molift Air or Molift Nomad lift motor and a sling, Molift Rail System Duo solves all overhead lifting situations, whatever the environment. With a lifting capacity of 300 kg (660 lbs), it is ideal for lifting and transferring a person from or to; bed, floor, chair, wheelchair, toilet, etc.

Suitable for all bed and room sizes

This 2 post gantry system has a rail length between 200 cm (78,7") and 350 cm (137,8") and has 4 steps adjustable height up to 255 cm (100,4"), making this freestanding hoist unit suitable for all bed and room sizes.

SPECIFICATIONS

Available rail lengths

200 - 250 - 300 - 350 cm

Height

225-255 cm

Depth	120 cm
Maximum internal space	350 x 240 cm
Safe Working Load SWL	300 kg
Material	Anodized aluminum, plastic
Total Weight	38-44,3 kg
Transport height	170 cm
Weight	Post, complete: 14.8 kg
	Boom: 4.2 kg/m
	Foot: 3.8 kg
	Post, without foot: 11 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au