




MOLIFT EVOSLING HYGIENE PADDED



SKU: 17222-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Molift EvoSling Hygiene Padded is made for people with strong core muscles and stability. It's designed to make daily tasks easier, safer, and more comfortable. Every part of this sling is made with the user in mind. It's not just about safety and comfort, but also about making everyday care routines simpler.

FEATURES

- **Large Leg-Opening:** Unique design with a large leg opening ensures easy access for toileting needs, removing the hassle from personal hygiene routines.
- **Ease of Use:** This sling is simple to apply and effortless to remove, saving time and ensuring user comfort.
- **Added Security:** Comes with a belt that guarantees additional security during lifting, ensuring the user feels safe and supported.
- **Optimal Positioning:** Designed to be user-friendly, it's a breeze to position, thanks to the application pockets on the leg supports.
- **Tailored Design:** The Molift EvoSling Hygiene Padded features a low back and large leg support opening, specifically crafted to facilitate dressing and undressing during toilet use.
- **Quality Material & Comfort:** Made with high-quality polyester, this sling boasts padded back and leg supports, ensuring users experience ultimate comfort.

SPECIFICATIONS

Item Variant	Item number	Height (mm)	Recommended user weight (kg)	Waist size (mm)
--------------	-------------	-------------	------------------------------	-----------------

XS	1722210	290	17-25	700-940
S	1722220	330	25-50	760-1000
M	1722230	330	45-95	850-1100
L	1722240	330	90-165	940-1200
XL	1722250	340	160-240	1020-1270
XXL	1722260	350	230-300	1090-1350

BROCHURE

[Molift Evosling Brochure](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
 1300 633 468

support@medimarnt.com.au