




## MOLIFT EVOSLING HYGIENE PADDED

**\$334.00 - \$413.00** Price range: \$334.00 through \$413.00



**SKU:** 17222-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Molift EvoSling Hygiene Padded is made for people with strong core muscles and stability. It's designed to make daily tasks easier, safer, and more comfortable. Every part of this sling is made with the user in mind. It's not just about safety and comfort, but also about making everyday care routines simpler.

### FEATURES

- **Large Leg-Opening:** Unique design with a large leg opening ensures easy access for toileting needs, removing the hassle from personal hygiene routines.
- **Ease of Use:** This sling is simple to apply and effortless to remove, saving time and ensuring user comfort.
- **Added Security:** Comes with a belt that guarantees additional security during lifting, ensuring the user feels safe and supported.
- **Optimal Positioning:** Designed to be user-friendly, it's a breeze to position, thanks to the application pockets on the leg supports.
- **Tailored Design:** The Molift EvoSling Hygiene Padded features a low back and large leg support opening, specifically crafted to facilitate dressing and undressing during toilet use.
- **Quality Material & Comfort:** Made with high-quality polyester, this sling boasts padded back and leg supports, ensuring users experience ultimate comfort.

## SPECIFICATIONS

Item Variant	Item number	Height (mm)	Recommended user weight (kg)	Waist size (mm)
--------------	-------------	-------------	------------------------------	-----------------

XS	1722210	290	17-25	700-940
S	1722220	330	25-50	760-1000
M	1722230	330	45-95	850-1100
L	1722240	330	90-165	940-1200
XL	1722250	340	160-240	1020-1270
XXL	1722260	350	230-300	1090-1350

## BROCHURE

[Molift Evosling Brochure](#)