



## **MOLIFT HARD TRAVEL CASE FOR SMART LIFTER 150**

**\$2,060.00**



**SKU: 990200**

### **PRODUCT DESCRIPTION**

### **PURPOSE AND BENEFITS**

The Molift Smart Travel Case makes travelling with the [Molift Smart 150](#) convenient and easy. The Molift Smart Travel Case has been designed uniquely with reinforced plastics, making it robust and protective against damage whilst your patient lifter is in transit. If you're travelling by plane, train or car this Travel Case will ensure that your Molift Smart 150, charger, battery and slings will arrive at their destination securely. For easy transportation the Molift Smart Hard Travel Case comes with wheels, allowing you to roll your Molift Smart 150 to any location.

To fit your Molift Smart 150 into the Travel Case you will need to fold it according to the operator manual, ensuring that the suspension is in the hook and the cord is around the castors and legs. Fit the battery and charger towards the bottom section.

For added security this Travel Case comes with a lock and key, TSA Approved for plane flights.

### **FEATURES**

- Designed uniquely with reinforced plastics
- Suitable for travel on airplanes
- TSA approved lock and key for air travel
- Fitted with wheels for easy transportation

### **SPECIFICATION**

**Colour:** Black

**Product Weight:** 11kg

**Size:**

Length 51" / 130cm

Width 22" / 55cm

Height 16" / 40cm