

MOLIFT MOVER 180 PATIENT LIFTER



SKU: 18000-2PT-COMBO-H



PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Molift Mover 180 is an all-round patient lifter suited to home care, aged care, and hospital use.

Molift Mover 180 patient lifter has a lightweight aluminum construction that makes this lifter easy to maneuver and easier on caregivers whilst offering a safe working load of 180kgs. A large lifting range, 27 – 163 cm (10.6" – 64.2") allows lifting from the highest beds as well as from the floor. This mobile patient lifter can also be supplied with 2 or 4 point sling bar.

Ergonomically correct: This patient lift is designed with both the user and the caregiver in mind. The ergonomically designed unique push bar handle provides flexible grip and multiple positions. The Mover 180 is powered by Molift's unique NiMH environmentally friendly batteries. This mobile patient lifter, like all Molift lifters, is equipped with a service notification system that lets you know when a service is due based on actual usage. When used in conjunction with the Molift service tool, data relating to lifter usage and patterns can be extracted.

FEATURES

- **Unique handle with many grip possibilities** - The unique handle frame provides the best natural grip, creating good ergonomics.
- **Integrated service software** - Automatic service notification system determines when the lifter requires service and maintenance.
- **Easy to maneuver** - The Molift Mover 180 has a total weight of only 32 kg / 70.5 lbs which, together with the lift's design, makes it easy to manoeuvre.
- **Easy lifting from chair** - Molift Mover 180 lifts easily from a chair, wheelchair, toilet and toilet chair.

- **Easy lift from a bed** - Molift Mover 180 lifts easily from bed with Molift Sling solution or the Molift Stretcher.

SPECIFICATIONS

Item number	M18000
Battery	NiMH 14,4 V-2,6 Ah
Battery charge time (hours)	3
Hoisting height max (mm)	1630
Hoisting height min (mm)	270
Hoisting interval (mm)	1365
Hoisting speed 75 kg (mm/sec)	50
Leg height (mm)	105
Length (mm)	1260
Material	Aluminium, steel
Maximum user weight (kg)	180
Rated performance	50 hoists at 75 kg
Service software	Molift Service Tool
Weight (kg)	32
Width chassis, legs open inner (mm)	1025
Width chassis, outer (mm)	670

SAFETY

- Automatic service notification system for lifter maintenance in Integrated service software
- Unique push bar for natural grip, good ergonomics, and various grip possibilities
- Compatible with 2-point and 4-point sling bars, with easily replaceable variants
- Lightweight and easy to maneuver

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA
 5048
 08 7079 8773
brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave,
 Elizabeth Vale SA 5112
 08 7079 8777
elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA
 5161
 08 7079 8771
reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au