



## MOLIFT MOVER 180 PATIENT LIFTER

**\$4,550.00 - \$4,600.00** Price range: \$4,550.00 through \$4,600.00



**SKU:** 18000-2PT-COMBO-H

## PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

The Molift Mover 180 is an all-round patient lifter suited to home care, aged care, and hospital use.

Molift Mover 180 patient lifter has a lightweight aluminum construction that makes this lifter easy to maneuver and easier on caregivers whilst offering a safe working load of 180kgs. A large lifting range, 27 - 163 cm (10.6" - 64.2") allows lifting from the highest beds as well as from the floor. This mobile patient lifter can also be supplied with 2 or 4 point sling bar.

**Ergonomically correct:** This patient lift is designed with both the user and the caregiver in mind. The ergonomically designed unique push bar handle provides flexible grip and multiple positions. The Mover 180 is powered by Molift's unique NiMH environmentally friendly batteries. This mobile patient lifter, like all Molift lifters, is equipped with a service notification system that lets you know when a service is due based on actual usage. When used in conjunction with the Molift service tool, data relating to lifter usage and patterns can be extracted.

### FEATURES

- **Unique handle with many grip possibilities** - The unique handle frame provides the best natural grip, creating good ergonomics.
- **Integrated service software** - Automatic service notification system determines when the lifter requires service and maintenance.
- **Easy to maneuver** - The Molift Mover 180 has a total weight of only 32 kg / 70.5 lbs which, together with the lift's design, makes it easy to manoeuver.
- **Easy lifting from chair** - Molift Mover 180 lifts easily from a chair, wheelchair, toilet and toilet chair.

- **Easy lift from a bed** - Molift Mover 180 lifts easily from bed with Molift Sling solution or the Molift Stretcher.

## SPECIFICATIONS

Item number	M18000
Battery	NiMH 14,4 V-2,6 Ah
Battery charge time (hours)	3
Hoisting height max (mm)	1630
Hoisting height min (mm)	270
Hoisting interval (mm)	1365
Hoisting speed 75 kg (mm/sec)	50
Leg height (mm)	105
Length (mm)	1260
Material	Aluminium, steel
Maximum user weight (kg)	180
Rated performance	50 hoists at 75 kg
Service software	Molift Service Tool
Weight (kg)	32
Width chassis, legs open inner (mm)	1025
Width chassis, outer (mm)	670

## SAFETY

- Automatic service notification system for lifter maintenance in Integrated service software
- Unique push bar for natural grip, good ergonomics, and various grip possibilities
- Compatible with 2-point and 4-point sling bars, with easily replaceable variants
- Lightweight and easy to maneuver