




## MOLIFT QUATTRO GANTRY FRAME



**SKU:** MQGF

## PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

#### A freestanding traverse gantry

Molift Quattro is a freestanding hoist unit that offers a (four post) Traverse System which is the best alternative when the hoisting need is temporary or when the facility is not build for fixed installations. The gantry is easy to install without the need to mount additional fixtures on wall or ceiling. It is ideally suited for temporary hoist applications in nursing or private homes.

#### Solves most overhead lifting situations

The traverse has a lightweight aluminum frame and four steel posts that provides excellent stability and covers the surface of the hoisting requirement. Molift Quattro solves all overhead lifting situations together with the Molift Air or Molift Nomad hoist motor and a sling. With a lifting capacity of max 300 kg, it is ideal for lifting and transferring a person from or to bed, floor, chair, wheelchair, toilet etc.

### TWO OPTIONS

Molift Quattro can be delivered in two sizes: 3x3 meters and 4x4 meters. The length of the posts are fixed and the height is 2.46 m. Molift Air or Molift Nomad with trolleys are ordered separately.

### SPECIFICATION

Item	Depth (mm)	Rail length (mm)	Weight (kg)	Width (mm)
------	------------	------------------	-------------	------------

<b>Molift Quattro 3</b>	3160	3000	111	3160
<b>Molift Quattro 4</b>	4160	4000	132	4160

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
 Brighton SA 5048  
 08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
 Ave, Elizabeth Vale SA 5112  
 08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
 Reynella SA 5161  
 08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
 03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
 03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
 1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)