



## **MOLIFT QUICK RAISER 205 SIT TO STAND HOIST**



**SKU:** M29000-H

## **PRODUCT DESCRIPTION**

### **Purpose & Benefits of Molift Quick Raiser 205 Sit to Stand Hoist**

The **Molift Quick Raiser 205 Sit to Stand Hoist** is an ergonomically designed, electrically powered hoist developed to assist with safe and comfortable standing transfers for users with limited mobility. Ideal for healthcare environments and home care settings, the Quick Raiser 205 supports active users during sit-to-stand movements, minimising strain on both the user and caregiver.

Its unique angled mast and low base design allow for easy access around furniture and close positioning to the user, while the electrically adjustable leg base ensures adaptability in tight spaces. The hoist features a comfortable, padded knee support and intuitive handset controls, enhancing both user comfort and caregiver efficiency during daily transfers.

### **Key Features**

- Designed for active sit-to-stand transfers in clinical or home care settings
- Angled mast and low base for optimal positioning and manoeuvrability
- Electrically powered lifting and leg base adjustment
- Ergonomic handle and intuitive handset for caregiver control
- Comfortable padded knee support for user stability

- Compact footprint allows for use in tight spaces
- Compatible with Molift sling range for tailored support
- Supports safe and efficient daily mobility assistance

## Specifications

### Dimensions

- Height: 1140 mm
- Length Adjustable Leg Option: 1320 mm
- Length Fixed Leg Option: 1165 mm
- Width Chassis, Outer: 630 mm
- Width Chassis, Legs Inner: 220 mm
- Width Chassis, Legs Open Inner: 590 mm
- Product weight: 42.5 kg / 44.5 kg

### Technical Specifications

- Rated Performance: 100 hoists/75 kg
- Expected Lifetime: 30,000 cycles or 10 years
- Battery: SLA 24V 2.9 Ah
- Battery Charge Time (hours): 5
- Protection Class: IP24
- Service Software: Molift Service Tool
- Material: Aluminium, Plastic, Steel

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)