


## MOLIFT RAISER PRO HEEL STRAP



**SKU:** 83211

## PRODUCT DESCRIPTION

### MOLIFT RAISER PRO HEEL STRAP - PURPOSE & BENEFITS

Optimize your Molift Raiser Pro Sit-to-Stand Transfer Platform with the Molift Raiser Pro Heel Strap. A meticulously designed accessory, this heel strap ensures that the user's feet remain in the correct, most secure position throughout the sit-to-stand transfer process.

### FEATURES

- **Foot Positioning Support:** The Strap has been created to provide unmatched support, ensuring the user's feet remain in the optimal position during transfers.
- **Easy Mounting Mechanism:** Designed for easy attachment, the heel strap effortlessly mounts over the turning knobs located on the sides of the Raiser Pro.
- **Adjustable Length:** Catering to the varied needs of users, the Heel Strap is equipped with two holes on each side, allowing for adjustment into three different lengths, ensuring a snug and comfortable fit for everyone.
- **Quality Material:** Crafted from soft silicone, the strap offers both durability and comfort, ensuring the user's heels feel supported without any discomfort.
- **Maintenance Simplicity:** With a design that emphasizes ease of use, this heel strap is not only simple to operate but also a breeze to clean.

**Please Note:** It's essential to lock the wheels of the Molift Raiser Pro before mounting the heel strap onto the turning knobs.

Make every sit-to-stand transfer with the Molift Raiser Pro a safer and more comfortable experience by incorporating the Molift Heel Strap. Prioritizing the user's safety and comfort, it's an indispensable addition to your transfer equipment.

## SPECIFICATIONS

- **Max User Weight:** 170kg
- **Material:** Silicone

You can purchase the Molift Raiser Pro - Sit-To-Stand aid with Medimart - [Click Here](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)