

## MOLIFT RAISER PRO HEEL STRAP

\$200.00



SKU: 83211



## PRODUCT DESCRIPTION

### MOLIFT RAISER PRO HEEL STRAP - PURPOSE & BENEFITS

Optimize your Molift Raiser Pro Sit-to-Stand Transfer Platform with the Molift Raiser Pro Heel Strap. A meticulously designed accessory, this heel strap ensures that the user's feet remain in the correct, most secure position throughout the sit-to-stand transfer process.

### FEATURES

- **Foot Positioning Support:** The Strap has been created to provide unmatched support, ensuring the user's feet remain in the optimal position during transfers.
- **Easy Mounting Mechanism:** Designed for easy attachment, the heel strap effortlessly mounts over the turning knobs located on the sides of the Raiser Pro.
- **Adjustable Length:** Catering to the varied needs of users, the Heel Strap is equipped with two holes on each side, allowing for adjustment into three different lengths, ensuring a snug and comfortable fit for everyone.
- **Quality Material:** Crafted from soft silicone, the strap offers both durability and comfort, ensuring the user's heels feel supported without any discomfort.
- **Maintenance Simplicity:** With a design that emphasizes ease of use, this heel strap is not only simple to operate but also a breeze to clean.

**Please Note:** It's essential to lock the wheels of the Molift Raiser Pro before mounting the heel strap onto the turning knobs.

Make every sit-to-stand transfer with the Molift Raiser Pro a safer and more comfortable experience by incorporating the Molift Heel Strap. Prioritizing the user's safety and comfort, it's an indispensable addition to your transfer equipment.

## SPECIFICATIONS

- **Max User Weight:** 170kg
- **Material:** Silicone

You can purchase the Molift Raiser Pro - Sit-To-Stand aid with Medimart - [Click Here](#)