



MOLIFT RAISER PRO HEEL STRAP

\$200.00 Original price was: \$200.00.\$180.00 Current price is: \$180.00.



SKU: 83211

PRODUCT DESCRIPTION

MOLIFT RAISER PRO HEEL STRAP - PURPOSE & BENEFITS

Optimize your Molift Raiser Pro Sit-to-Stand Transfer Platform with the Molift Raiser Pro Heel Strap. A meticulously designed accessory, this heel strap ensures that the user's feet remain in the correct, most secure position throughout the sit-to-stand transfer process.

FEATURES

- **Foot Positioning Support**: The Strap has been created to provide unmatched support, ensuring the user's feet remain in the optimal position during transfers.
- **Easy Mounting Mechanism**: Designed for easy attachment, the heel strap effortlessly mounts over the turning knobs located on the sides of the Raiser Pro.
- **Adjustable Length**: Catering to the varied needs of users, the Heel Strap is equipped with two holes on each side, allowing for adjustment into three different lengths, ensuring a snug and comfortable fit for everyone.
- **Quality Material**: Crafted from soft silicone, the strap offers both durability and comfort, ensuring the user's heels feel supported without any discomfort.
- **Maintenance Simplicity**: With a design that emphasizes ease of use, this heel strap is not only simple to operate but also a breeze to clean.

Please Note: It's essential to lock the wheels of the Molift Raiser Pro before mounting the heel strap onto the turning knobs.

Make every sit-to-stand transfer with the Molift Raiser Pro a safer and more comfortable experience by incorporating the Molift Heel Strap. Prioritizing the user's safety and comfort, it's an indispensable addition to your transfer equipment.



SPECIFICATIONS

• Max User Weight: 170kg

• Material: Silicone

You can purchase the Molift Raiser Pro - Sit-To-Stand aid with Medimart - Click Here