


MOLIFT RAISER PRO STRAP WITH SLIDING SLEEVE

\$385.00



SKU: 82528-H



PRODUCT DESCRIPTION

MOLIFT RAISER PRO STRAP WITH SLIDING SLEEVE - PURPOSE & BENEFITS

Upgrade your Molift Raiser Pro Sit to Stand Transfer Platform with the Molift Raiser Pro Strap with Sliding Sleeve. Expertly crafted to offer both security and convenience, this safety strap ensures that sit-to-stand transfers are smoother and more comfortable for the user.

FEATURES

- **Enhanced Support:** For those who need extra assistance during sit-to-stand transfers, this strap provides much-needed support, ensuring each transfer is smooth and secure.
- **Innovative Sliding Design:** Whether you have one or two caregivers assisting, the strap is designed to safely slide, letting caregivers effortlessly pull the user upward from a seated position.
- **Secure Fastening:** One end of the strap attaches seamlessly to the designated attachment points on the side of the handle. The other end is skillfully pulled by the caregiver, ensuring the user is gently glided to a standing posture.
- **Size Variability:** To cater to a broader range of users, the strap is available in two distinct sizes: Small/Medium and Large/Extra Large. Please make your desired selection from the dropdown menu above.
- **User-Friendly & Hygienic:** The strap is not only easy to use but also easy to maintain. It's machine washable and can be tumble dried at a gentle temperature. For those keen on a deeper cleanse, the strap can handle temperatures up to 80°C.

With the Molift Raiser Pro Strap With Sliding Sleeve, add an extra layer of safety, support, and ease to your transfer routines. Prioritize comfort and security with this essential accessory.

SPECIFICATIONS

- **Sliding Sleeve Length:** S/M = 79cm, L/XL = 107cm
- **Length:** S/M = 153cm, L/XL = 193cm

You can purchase the Molift Raiser Pro - Sit-To-Stand aid with Medimart - [Click Here](#)