


## MOLIFT RAISER PRO STRAP WITH SLIDING SLEEVE



**SKU:** 82528-H



## PRODUCT DESCRIPTION

### MOLIFT RAISER PRO STRAP WITH SLIDING SLEEVE - PURPOSE & BENEFITS

Upgrade your Molift Raiser Pro Sit to Stand Transfer Platform with the Molift Raiser Pro Strap with Sliding Sleeve. Expertly crafted to offer both security and convenience, this safety strap ensures that sit-to-stand transfers are smoother and more comfortable for the user.

## FEATURES

- **Enhanced Support:** For those who need extra assistance during sit-to-stand transfers, this strap provides much-needed support, ensuring each transfer is smooth and secure.
- **Innovative Sliding Design:** Whether you have one or two caregivers assisting, the strap is designed to safely slide, letting caregivers effortlessly pull the user upward from a seated position.
- **Secure Fastening:** One end of the strap attaches seamlessly to the designated attachment points on the side of the handle. The other end is skillfully pulled by the caregiver, ensuring the user is gently glided to a standing posture.
- **Size Variability:** To cater to a broader range of users, the strap is available in two distinct sizes: Small/Medium and Large/Extra Large. Please make your desired selection from the dropdown menu above.
- **User-Friendly & Hygienic:** The strap is not only easy to use but also easy to maintain. It's machine washable and can be tumble dried at a gentle temperature. For those keen on a deeper cleanse, the strap can handle temperatures up to 80°C.

With the Molift Raiser Pro Strap With Sliding Sleeve, add an extra layer of safety, support, and ease to your transfer routines. Prioritize comfort and security with this essential accessory.

## SPECIFICATIONS

- **Sliding Sleeve Length:** S/M = 79cm, L/XL = 107cm
- **Length:** S/M = 153cm, L/XL = 193cm

You can purchase the Molift Raiser Pro - Sit-To-Stand aid with Medimart - [Click Here](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)