



## **MOLIFT RAISER PRO - SIT TO STAND TRANSFER PLATFORM**

**\$2,009.00**



**SKU: 16090204**

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

The Molift Raiser Pro is an ergonomic platform for a safe and active sit-to-stand transfer. The design imitates the natural pattern of movement and fulfils the transfer needs of users with some weight bearing capacity. The handle has two different fixing heights for accessories and the soft padded leg supports can be adjusted.

### **FEATURES**

- Ergonomics for the user
- Natural Movement Patterns
- Transfers in Confined Spaces
- Works in all Healthcare Environments
- Versatile Handle
- Central Brake
- Natural Movement Patterns
- Safe and Efficient
- Low Step-On Height
- Optional accessories: Soft handles, Raiser strap, Raiser belt

## **SPECIFICATIONS**

- Height: 1200 mm
- Depth: 630 mm
- Width: 54 mm
- Step-on Height: 35 mm
- Total Weight: 10.7 kg
- Weight of Baseplate: 5.4 kg
- Maximum User Weight: 170 kg
- Material: Powder-coated aluminium, plastic, steel, chrome

## **ATTACHMENTS**

[Click here to view the attachments](#)