



MOLIFT RAISER PRO - SIT TO STAND TRANSFER PLATFORM

\$2,009.00 Original price was: \$2,009.00.<u>\$1,808.00</u>Current price is: \$1,808.00.



SKU: 16090204

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Molift Raiser Pro is an ergonomic platform for a safe and active sit-to-stand transfer. The design imitates the natural pattern of movement and fulfils the transfer needs of users with some weight bearing capacity. The handle has two different fixing heights for accessories and the soft padded leg supports can be adjusted.

FEATURES

- Ergonomics for the user
- Natural Movement Patterns
- Transfers in Confined Spaces
- Works in all Healthcare Environments
- Versatile Handle
- Central Brake
- Natural Movement Patterns
- Safe and Efficient
- Low Step-On Height
- Optional accessories: Soft handles, Raiser strap, Raiser belt



SPECIFICATIONS

Height: 1200 mmDepth: 630 mmWidth: 54 mm

Step-on Height: 35 mmTotal Weight: 10.7 kg

Weight of Baseplate: 5.4 kgMaximum User Weight: 170 kg

• Material: Powder-coated aluminium, plastic, steel, chrome

ATTACHMENTS

Click here to view the attachments