




MOLIFT RGOSLING ACTIVE

\$394.00 - \$620.00 Price range: \$394.00 through \$620.00



SKU: 17209-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Molift RgoSling Active, the embodiment of a sophisticated and efficient hoisting solution specifically designed for active transfers using the Molift Quick Raiser 205.

FEATURES

- **Specially Designed:** Crafted for active hoisting, this sling seamlessly pairs with the Molift Quick Raiser 205. It supports the lumbar region perfectly, making transitions between different seated positions a breeze.
- **Comfortable Padded Support:** Experience unparalleled comfort with its padded structure. This ensures not just safety but a snug fit for the user during all transfers.
- **Self-adjusting Sliding Loops:** The RgoSling Active comes equipped with sliding loops that adjust on their own to fit the user's size, ensuring a secure grip and added support around the hips during stand-up procedures.
- **4-Point Active Sling Bar:** This feature ensures a balanced and stable hoisting experience, catering to a range of user needs.
- **Effortless Transfers:** Designed for a hassle-free experience, the RgoSling Active is simple to place around the user. Its efficient anti-slip lining ensures the sling remains in place, preventing upward movement during transfers. Furthermore, its wide back support evenly distributes the user's weight, maximizing comfort. Additionally, the waist belt is an added layer of support, providing stability during movements.
- **Active Hoisting:** Elevate your hoisting experience with this state-of-the-art feature, enhancing the comfort and safety of every transfer.

Step into the future of active transfers with the Molift RgoSling Active, where every design element

focuses on user comfort, safety, and ease of use. Experience the revolution in transfer support today!

SPECIFICATIONS

Item Variant	Item number	Height (mm)	Waist size (mm)
XS	1720910	370	500 - 700
S	1720920	410	600 - 800
M	1720930	420	750 - 950
L	1720940	420	950 - 1200
XL	1720950	420	1200 - 1400

Buy the Molift RgoSling Active with the [NDIS payment method](#)