



MOLIFT RGOSLING AMBULATING VEST

\$460.00 - \$600.00 Price range: \$460.00 through \$600.00



SKU: 17213-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Crafted with precision and inspired by the best of climbing technology, the Molift RgoSling Ambulating Vest offers an unparalleled experience for those in rehabilitation. From gait training to balance exercises, this vest is your ideal companion.

FEATURES

- **Multifunctional Use:** Whether you're embarking on gait training, standing exercises, or honing your balance, this vest supports both independent and assisted training.
- **Superior Comfort & Fit:** With its double-buckled, wide padded belt, users can be assured of even pressure distribution across a larger area, ensuring maximum comfort throughout the rehabilitation process.
- **Secure Positioning:** The optional padded groin straps are ingeniously designed to prevent the vest from sliding up, keeping it securely in place during use.
- **Climbing Inspired Design:** Borrowing key features from professional climbing gear, the Molift RgoSling Ambulating Vest prioritizes user comfort and safety. The self-tightening suspension loops serve as a protective measure against slips or falls.
- **Independent Training Boost:** The front buckle system broadens the horizon for users, empowering them with the potential to engage in gait training or toileting procedures, possibly even without any assistance.

SPECIFICATIONS

Item Variant	Item number	Height (mm)	Recommended user weight (kg)	Waist size (mm)
--------------	-------------	-------------	------------------------------	-----------------

XXS	1721300	150	12 - 17	500 - 650
XS	1721310	170	17 - 25	600 - 750
S	1721320	190	25 - 50	700 - 900
M	1721330	210	45 - 95	850 - 1050
L	1721340	220	90 - 165	1000 - 1300
XL	1721350	230	160 - 240	1250 - 1550