



## MOLIFT RGOSLING AMPUTEE HIGHBACK PADDED



**SKU:** 17212-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Designed with precision and tailored to meet the unique needs of users with an amputated leg or those with diminished muscle tone, the Molift RgoSling Amputee HighBack Padded offers enhanced safety, comfort, and functionality during hoisting and transferring situations.

### FEATURES

- **Specialized Support:** Created especially for users with amputated limbs or compromised muscle tone, ensuring that every lift or transfer is secure and comfortable.
- **HighBack Design:** For those needing added trunk and head support, the HighBack feature offers a snug fit, ensuring maximum stability and confidence during every transfer.
- **Safety Waist Belt:** With a built-in waist belt, safety takes the forefront, giving an extra layer of security during lifts.
- **Narrow Leg Support:** The sling's design includes a narrow leg support opening that wraps securely under both legs, minimizing any risk during hoisting. This feature provides heightened safety by preventing any accidental slips or movements.
- **Extended Leg Application:** The divided leg parts come with extra length, allowing them to be easily positioned under both legs. With additional straps integrated into the design, there's no risk of the legs slipping out of place, further boosting user confidence during transfers.
- **Anatomical Cut & Padding:** Every curve and padding on the RgoSling Amputee HighBack has been meticulously designed to provide unbeatable support throughout the entirety of the hoisting and transferring process.
- **Versatile Use:** Whether you're moving from a horizontal position or adjusting from sitting, the Molift RgoSling Amputee HighBack Padded caters to various hoisting needs, making it a

versatile choice for caregivers and users alike.

## SPECIFICATIONS

Item Variant	Item number	Height (mm)	Material	Recommended user weight (kg)	Waist size (mm)
XXS	1721200	570	Polyester, PE-cellfoam	12-17	600-850
XS	1721210	690	Polyester, PE-cellfoam	17-25	670-920
S	1721220	790	Polyester, PE-cellfoam	25-50	750-1050
M	1721230	920	Polyester, PE-cellfoam	45-95	900-1100
L	1721240	980	Polyester,, PE-cellfoam	90-165	1100-1350
XL	1721250	1000	Polyester, PE-cellfoam	160-240	1120-1400
XXL	1721260	1000	Polyester, PE-cellfoam	230-300	1150-1500

Buy the Molift RgoSling Amputee Highback Padded with the [NDIS payment method](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)