



MOLIFT RGOSLING AMPUTEE MEDIUMBACK PADDED



SKU: 17211-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Tailored specifically for users with an amputated leg or those with reduced muscle tone, the Molift RgoSling Amputee Mediumback Padded is the embodiment of safety and comfort. Designed with keen attention to detail, this sling offers an unparalleled hoisting and transferring experience.

FEATURES

- **Specialized Support:** Crafted with the unique needs of users with amputated limbs or diminished muscle tone in mind, ensuring optimal comfort and safety during each transfer.
- **MediumBack Design:** Balancing support and flexibility, the MediumBack design ensures users have the necessary trunk support without compromising movement.
- **Safety First with Waist Belt:** The built-in waist belt offers an added layer of security, enhancing the user's confidence during lifts.
- **Secure Leg Support:** Featuring a narrow leg support opening, the design snugly wraps under both legs, ensuring minimal movement and maximum safety during hoists.
- **Extended Leg Application:** The leg parts are not only divided for ease but also come with extra length, ensuring they fit comfortably under both legs. The added straps further ensure the legs stay in position, enhancing the user's trust in the product.
- **Anatomical Cut & Padding:** Every contour and padding on the Molift RgoSling Amputee Mediumback Padded has been meticulously designed. This attention to detail ensures unbeatable support throughout the hoisting and transferring process, making every movement seamless.

PURPOSE & BENEFITS

Item Variant	Item number	Height (mm)	Recommended user weight (kg)	Waist size (mm)	Waist size (mm)
XXS	1721100	440	12 - 17	600-850	600-850
XS	1721110	510	17 - 25	670-920	670-920
S	1721120	570	25-50	750-1050	750-1050
M	1721130	620	45 - 95	900-1100	900-1100
L	1721140	640	90 - 165	1100-1350	1100-1350
XL	1721150	680	160 - 240	1120-1400	1120-1400
XXL	1721160	680	230 - 300	1150-1500	1150-1500

Buy the Molift RgoSling Amputee Mediumback Padded with the [NDIS payment method](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au

