



MOLIFT RGOSLING FABRIC STRETCHER



SKU: 17215-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Molift RgoSling Fabric Stretcher, a soft fabric hoisting device meticulously designed for horizontal and lateral transfers. This innovative device is ideal for hospitals and institutions where users require a stable recumbent position due to medical conditions.

FEATURES

- **Soft Fabric Hoisting:** Crafted with gentle fabric, the RgoSling Fabric Stretcher makes hoisting a comfortable and smooth process, bringing a touch of care to the routine.
- **Versatile Transfers:** Whether you're dealing with horizontal or lateral transfers, the Fabric Stretcher caters to a range of needs, making it a versatile choice in medical environments.
- **Stretcher Alternative:** For those situations where using a traditional stretcher is inappropriate, such as for comatose or sedated clients, the Fabric Stretcher provides a comfortable and safer alternative.
- **8-Point Sling Bar:** This product necessitates an 8-point sling bar for hoisting. This bar is compatible with Molift Partner 255 and Molift Air.
- **Integrated Hoisting Straps:** The Molift Fabric Stretcher comes with integrated hoisting straps that hook directly to the sling bar. The straps offer some degree of adjustability in length through the use of loops. This allows for an open and stable horizontal hoist, enhancing the safety and comfort of the procedure.

SPECIFICATIONS

Item Variant	Item number	Length (mm)	Width (mm)
--------------	-------------	-------------	------------

M	1721520	1750	450
L	1721530	1780	650
XL	1721540	1830	900

Buy the Molift RgoSling Fabric Stretcher with the [NDIS payment method](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au