



## MOLIFT RGOSLING FABRIC STRETCHER

**\$780.00 - \$1,080.00** Price range: \$780.00 through \$1,080.00



**SKU:** 17215-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Introducing the Molift RgoSling Fabric Stretcher, a soft fabric hoisting device meticulously designed for horizontal and lateral transfers. This innovative device is ideal for hospitals and institutions where users require a stable recumbent position due to medical conditions.

### FEATURES

- **Soft Fabric Hoisting:** Crafted with gentle fabric, the RgoSling Fabric Stretcher makes hoisting a comfortable and smooth process, bringing a touch of care to the routine.
- **Versatile Transfers:** Whether you're dealing with horizontal or lateral transfers, the Fabric Stretcher caters to a range of needs, making it a versatile choice in medical environments.
- **Stretcher Alternative:** For those situations where using a traditional stretcher is inappropriate, such as for comatose or sedated clients, the Fabric Stretcher provides a comfortable and safer alternative.
- **8-Point Sling Bar:** This product necessitates an 8-point sling bar for hoisting. This bar is compatible with Molift Partner 255 and Molift Air.
- **Integrated Hoisting Straps:** The Molift Fabric Stretcher comes with integrated hoisting straps that hook directly to the sling bar. The straps offer some degree of adjustability in length through the use of loops. This allows for an open and stable horizontal hoist, enhancing the safety and comfort of the procedure.

### SPECIFICATIONS

Item Variant	Item number	Length (mm)	Width (mm)
--------------	-------------	-------------	------------

M	1721520	1750	450
L	1721530	1780	650
XL	1721540	1830	900

Buy the Molift RgoSling Fabric Stretcher with the [NDIS payment method](#)