



MOLIFT RGOSLING HIGHBACK PADDED

\$640.00 - \$780.00 Price range: \$640.00 through \$780.00



SKU: 17202-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Presenting the Molift RgoSling HighBack Padded, a versatile all-around sling built with ergonomic considerations for comfortable hoisting. Perfectly suited for most users and lifting situations, it's an ideal choice for both homecare and institutional environments.

FEATURES

- **Ergonomic Comfort:** Designed with built-in ergonomics, this sling ensures a comfortable hoisting experience while providing comprehensive support for the entire body, including the head.
- **Flexible Usage:** Whether it's hoisting from a lying or sitting position, or even lifting from the floor, the RgoSling HighBack Padded is up to the task, making it a valuable addition to any care environment.
- **Support for the Head:** Not only does the HighBack provide body support, but it also supports the head. It allows the user to maintain a slightly inclined position, making it an excellent fit for individuals with compromised trunk and head stability.
- **Comfortable Padding:** With padded back and leg rests, this sling goes the extra mile to ensure comfort during each hoisting operation.
- **Versatile Application:** The Molift RgoSling Highback Padded fits most users and hoisting situations, making it a practical choice for a wide range of home care and institutional settings.

SPECIFICATIONS

Item Variant	Item number	Height (mm)	Leg support opening width (mm)	Maximum user weight (kg)	Recommended user weight (kg)
XXS	1720200	580	180	300	12-17
XS	1720210	680	200	300	17-25
S	1720220	800	270	300	25-50
M	1720230	920	320	300	45-95
L	1720240	980	400	300	90-160
XL	1720250	980	520	350	160-240
XXL	1720260	1000	660	350	230-350

Buy the Molift RgoSling Highback Padded with the [NDIS payment method](#)