



## MOLIFT RGOSLING MEDIUMBACK NET

**\$330.00 - \$550.00** Price range: \$330.00 through \$550.00



**SKU:** 17203-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Introducing the Molift RgoSling MediumBack Net, an ergonomic sling designed for comfortable and safe lifting during bathing and shower situations. This versatile bath sling is perfect for a variety of lifting situations in both homecare and institutional environments.

### FEATURES

- **Full Body Support:** The RgoSling MediumBack Net provides complete body support, ensuring user safety and comfort during lifting and transfers. It's especially suited for users who have good head control.
- **Breathable & Quick-Drying Material:** Crafted from soft net mesh polyester, this sling is not only durable and flexible but also boasts good ventilation. Its quick-drying feature makes it particularly handy in bathing and shower situations.
- **Ergonomic Design:** Designed with user comfort in mind, the MediumBack Net features an ergonomic design that ensures a safe and comfortable lift every time.
- **Similarities with Padded Versions:** The MediumBack Net retains the same shape and key features as its padded counterparts, including sliding loops, an application pocket in the back, and leg rest and straps for increased weight distribution in the leg parts.

### SPECIFICATIONS

Item Variant	Item number	Height (mm)	Leg support opening width (mm)	Recommended user weight (kg)
XXS	1720300	440	180	12-17

XS	1720310	500	200	17-25
S	1720320	560	270	25-50
M	1720330	620	320	45-95
L	1720340	650	400	90-165
XL	1720350	680	520	160-240
XXL	1720360	680	660	230-300

Buy the Molift RgoSling Mediumback Net with the [NDIS payment method](#)