



MOLIFT RGOSLING MEDIUMBACK NET



SKU: 17203-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Molift RgoSling MediumBack Net, an ergonomic sling designed for comfortable and safe lifting during bathing and shower situations. This versatile bath sling is perfect for a variety of lifting situations in both homecare and institutional environments.

FEATURES

- **Full Body Support:** The RgoSling MediumBack Net provides complete body support, ensuring user safety and comfort during lifting and transfers. It's especially suited for users who have good head control.
- **Breathable & Quick-Drying Material:** Crafted from soft net mesh polyester, this sling is not only durable and flexible but also boasts good ventilation. Its quick-drying feature makes it particularly handy in bathing and shower situations.
- **Ergonomic Design:** Designed with user comfort in mind, the MediumBack Net features an ergonomic design that ensures a safe and comfortable lift every time.
- **Similarities with Padded Versions:** The MediumBack Net retains the same shape and key features as its padded counterparts, including sliding loops, an application pocket in the back, and leg rest and straps for increased weight distribution in the leg parts.

SPECIFICATIONS

Item Variant	Item number	Height (mm)	Leg support opening width (mm)	Recommended user weight (kg)
XXS	1720300	440	180	12-17

XS	1720310	500	200	17-25
S	1720320	560	270	25-50
M	1720330	620	320	45-95
L	1720340	650	400	90-165
XL	1720350	680	520	160-240
XXL	1720360	680	660	230-300

Buy the Molift RgoSling Mediumback Net with the [NDIS payment method](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au