



MOLIFT RGOSLING MEDIUMBACK PADDED



SKU: 17201-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Welcome to the Molift RgoSling MediumBack Padded, an all-around sling designed for both comfort and utility. This sling is suitable for most users and a wide range of hoisting situations, making it a perfect fit for homecare and institutional environments.

FEATURES

- **Ergonomic Comfort:** The RgoSling MediumBack Padded is built with an ergonomic design that ensures comfortable hoisting. Its design supports the entire body, including users with impaired trunk and head stability.
- **Flexible Use:** This sling is a basic model designed for users who have good head control, making it adaptable for a variety of users and lifting situations.
- **Comfortable Padding:** It features a comfortable padded back and leg rests, ensuring a soft and secure experience during hoisting operations.
- **Versatile Application:** Whether used in home care or institutional environments, the RgoSling MediumBack Padded caters to most users and hoisting situations, making it a highly versatile choice for care providers.

SPECIFICATIONS

Item Variant	Item number	Height (mm)	Leg support opening width (mm)	Maximum user weight (kg)
XXS	1720100	440	180	300
XS	1720110	500	200	300

S	1720120	560	270	300
M	1720130	620	320	300
L	1720140	650	400	300
XL	1720150	680	520	350
XXL	1720160	680	660	350

Buy the Molift RgoSling Mediumback Padded with the [NDIS payment method](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au