



## MOLIFT RGOSLING MEDIUMBACK PADDED

**\$512.00 - \$690.00** Price range: \$512.00 through \$690.00



**SKU:** 17201-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Welcome to the Molift RgoSling MediumBack Padded, an all-around sling designed for both comfort and utility. This sling is suitable for most users and a wide range of hoisting situations, making it a perfect fit for homecare and institutional environments.

### FEATURES

- **Ergonomic Comfort:** The RgoSling MediumBack Padded is built with an ergonomic design that ensures comfortable hoisting. Its design supports the entire body, including users with impaired trunk and head stability.
- **Flexible Use:** This sling is a basic model designed for users who have good head control, making it adaptable for a variety of users and lifting situations.
- **Comfortable Padding:** It features a comfortable padded back and leg rests, ensuring a soft and secure experience during hoisting operations.
- **Versatile Application:** Whether used in home care or institutional environments, the RgoSling MediumBack Padded caters to most users and hoisting situations, making it a highly versatile choice for care providers.

### SPECIFICATIONS

Item Variant	Item number	Height (mm)	Leg support opening width (mm)	Maximum user weight (kg)
XXS	1720100	440	180	300
XS	1720110	500	200	300

S	1720120	560	270	300
M	1720130	620	320	300
L	1720140	650	400	300
XL	1720150	680	520	350
XXL	1720160	680	660	350

Buy the Molift RgoSling Mediumback Padded with the [NDIS payment method](#)