




MOLIFT RGOSLING NET BATHING HIGHBACK

\$380.00 - \$680.00 Price range: \$380.00 through \$680.00



SKU: 17204-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Meet the Molift RgoSling Net Bathing Highback, an all-around, breathable sling designed for most users and hoisting situations, making it a versatile choice for both homecare and institutional settings.

FEATURES

- **Ergonomic Design:** Built with ergonomics in mind, this sling ensures comfortable hoisting while providing full-body and head support. It's perfect for users who need to transition from a lying to a sitting position and vice versa.
- **Head Support:** The HighBack design not only provides support to the body but also to the head. It facilitates a slightly inclined position, making it ideal for users with compromised trunk and head stability.
- **Breathable & Fast-Drying Fabric:** Made from soft net mesh polyester material, the RgoSling HighBack Net offers excellent ventilation, ensuring user comfort during use. Its fast-drying properties make it particularly suitable for shower and bathing situations.
- **Enhanced Weight Distribution:** Like its padded counterparts, the HighBack Net also features sliding loops, an application pocket in the back, and leg rests and straps for improved weight distribution in the leg parts.
- **Durable & Flexible:** Crafted with durability in mind, the flexible polyester net is designed to withstand regular use while maintaining its shape and functionality.

SPECIFICATIONS

Item Variant	Item number	Height (mm)	Leg support opening width (mm)	Recommended user weight (kg)
XXS	1720400	580	180	12-17
XS	1720410	680	200	17-25
S	1720420	800	270	25-50
M	1720430	920	320	45-95
L	1720440	980	400	90-165
XL	1720450	980	520	160-240
XXL	1720460	1000	660	230-300

Buy the Molift RgoSling Highback Net Suits Bathing with the [NDIS payment method](#)