




## MOLIFT RGOSLING NET BATHING HIGHBACK

**\$380.00 - \$680.00**



**SKU: 17204-H**

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Meet the Molift RgoSling Net Bathing Highback, an all-around, breathable sling designed for most users and hoisting situations, making it a versatile choice for both homecare and institutional settings.

### FEATURES

- **Ergonomic Design:** Built with ergonomics in mind, this sling ensures comfortable hoisting while providing full-body and head support. It's perfect for users who need to transition from a lying to a sitting position and vice versa.
- **Head Support:** The HighBack design not only provides support to the body but also to the head. It facilitates a slightly inclined position, making it ideal for users with compromised trunk and head stability.
- **Breathable & Fast-Drying Fabric:** Made from soft net mesh polyester material, the RgoSling HighBack Net offers excellent ventilation, ensuring user comfort during use. Its fast-drying properties make it particularly suitable for shower and bathing situations.
- **Enhanced Weight Distribution:** Like its padded counterparts, the HighBack Net also features sliding loops, an application pocket in the back, and leg rests and straps for improved weight distribution in the leg parts.
- **Durable & Flexible:** Crafted with durability in mind, the flexible polyester net is designed to withstand regular use while maintaining its shape and functionality.

## SPECIFICATIONS

Item Variant	Item number	Height (mm)	Leg support opening width (mm)	Recommended user weight (kg)
XXS	1720400	580	180	12-17
XS	1720410	680	200	17-25
S	1720420	800	270	25-50
M	1720430	920	320	45-95
L	1720440	980	400	90-165
XL	1720450	980	520	160-240
XXL	1720460	1000	660	230-300

Buy the Molift RgoSling Highback Net Suits Bathing with the [NDIS payment method](#)