




MOLIFT RGOSLING SHADOW HIGHBACK

\$380.00 - \$460.00 Price range: \$380.00 through \$460.00



SKU: 17229-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Molift RgoSling Shadow HighBack Net, a specially designed sling for users who require total support. It is uniquely constructed to be left under the user even after transfer, ensuring continuous comfort and ease of use. The absence of straps or seams on the user-contact surface makes it an optimal choice for individuals using moulded seat units.

FEATURES

- **Persistent Convenience:** The RgoSling Shadow HighBack Net is designed to stay with the user even when they are seated in their wheelchair. This feature allows for easy, ongoing transfers without the need for continuous sling removal and reapplication.
- **Full Body Support:** This sling offers comprehensive support for the entire core and head. Its high back and head support cater to users needing total support during transfers.
- **Versatile Lifting:** The RgoSling Shadow HighBack Net is suitable for various lifting scenarios, ranging from lying to sitting, or lying positions. This versatility ensures it can accommodate different user needs.
- **Suitable for Moulded Seats:** The sling is specially designed to be optimal for users with moulded seat units, offering a smooth, comfortable surface free of straps or seams.
- **Split Leg Support:** The sling features split leg support with a medium opening. However, it is not suitable for double-leg amputees.

SPECIFICATIONS

Item Variant	Item number	Height (mm)	Leg support opening width (mm)	Recommended user weight (kg)
--------------	-------------	-------------	--------------------------------	------------------------------

XXS	1722900	550	180	12-17
XS	1722910	680	200	17-25
S	1722920	800	270	25-50
M	1722930	920	320	45-95
L	1722940	980	400	90-165
XL	1722950	980	520	160-240