




MOLIFT RGOSLING SHADOW HIGHBACK



SKU: 17229-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Molift RgoSling Shadow HighBack Net, a specially designed sling for users who require total support. It is uniquely constructed to be left under the user even after transfer, ensuring continuous comfort and ease of use. The absence of straps or seams on the user-contact surface makes it an optimal choice for individuals using moulded seat units.

FEATURES

- **Persistent Convenience:** The RgoSling Shadow HighBack Net is designed to stay with the user even when they are seated in their wheelchair. This feature allows for easy, ongoing transfers without the need for continuous sling removal and reapplication.
- **Full Body Support:** This sling offers comprehensive support for the entire core and head. Its high back and head support cater to users needing total support during transfers.
- **Versatile Lifting:** The RgoSling Shadow HighBack Net is suitable for various lifting scenarios, ranging from lying to sitting, or lying positions. This versatility ensures it can accommodate different user needs.
- **Suitable for Moulded Seats:** The sling is specially designed to be optimal for users with moulded seat units, offering a smooth, comfortable surface free of straps or seams.
- **Split Leg Support:** The sling features split leg support with a medium opening. However, it is not suitable for double-leg amputees.

SPECIFICATIONS

Item Variant	Item number	Height (mm)	Leg support opening width (mm)	Recommended user weight (kg)
--------------	-------------	-------------	--------------------------------	------------------------------

XXS	1722900	550	180	12-17
XS	1722910	680	200	17-25
S	1722920	800	270	25-50
M	1722930	920	320	45-95
L	1722940	980	400	90-165
XL	1722950	980	520	160-240

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
 1300 633 468

support@medimarnt.com.au