



MOLIFT RGOSLING STANDUP PADDED



SKU: 17207-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Molift RgoSling StandUp Padded, the perfect sling solution for active transfers with Molift Quick Raiser. This product is specifically designed to offer comfort and lumbar support during hoisting and transition between different sitting positions.

FEATURES

- **Optimized for Active Hoisting:** This sling is specifically designed for use with the Molift Quick Raiser, ensuring seamless and comfortable transfers.
- **Ideal for Toileting Situations:** With its easy placement and secure design, the RgoSling StandUp Padded is an excellent choice for use in toileting situations.
- **Compatibility with 2-Point Sling Bar:** Designed for use with a 2-point sling bar, this product is easy to install and use.
- **Enhanced Comfort and Stability:** The sling features an efficient anti-slip lining to prevent sliding during use. Its wide back support evenly distributes weight and pressure over a large surface, enhancing user comfort. An additional waist belt provides extra support and stability.
- **Additional Support Option:** For users requiring extra pelvic support, the Molift RgoSling StandUp Padded with support variant is an ideal alternative.

SPECIFICATIONS

Item Variant	Item number	Height (mm)	Waist size (mm)
XS	1720710	300	550-700
S	1720720	340	650-850

M	1720730	350	750-950
L	1720740	350	900-1200
XL	1720750	350	1150-1450

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
 1300 633 468

support@medimarnt.com.au