



## MOLIFT RGOSLING STANDUP PADDED WITH SUPPORT

**\$649.00 - \$786.00** Price range: \$649.00 through \$786.00



**SKU:** 17208-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Introducing the Molift RgoSling StandUp Padded With Support, the perfect sling solution for active transfers with Molift Quick Raiser. This product is specifically designed to offer comfort and lumbar support during hoisting and transition between different sitting positions.

### FEATURES

- **Optimized for Active Hoisting:** This sling is specifically designed for use with the Molift Quick Raiser, ensuring seamless and comfortable transfers.
- **Ideal for Toileting Situations:** With its easy placement and secure design, the RgoSling StandUp Padded is an excellent choice for use in toileting situations.
- **Compatibility with 2-Point Sling Bar:** Designed for use with a 2-point sling bar, this product is easy to install and use.
- **Enhanced Comfort and Stability:** The sling features an efficient anti-slip lining to prevent sliding during use. Its wide back support evenly distributes weight and pressure over a large surface, enhancing user comfort. An additional waist belt provides extra support and stability.

### SPECIFICATIONS

Item Variant	Item number	Height (mm)	Waist size (mm)
XS with Support	1720810	300	550-700
S with Support	1720820	340	650-850
M with Support	1720830	350	750-950
L with Support	1720840	350	900-1200

XL with Support	1720850	350	1150-1450
-----------------	---------	-----	-----------