



MOLIFT RGOSLING TOILET HIGHBACK PADDED



SKU: 17206-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Molift RgoSling Toilet HighBack, a highback flexible sling engineered to provide optimal functions for toilet situations. This product is designed with large leg openings, making it easy to apply, remove, and facilitate toileting. With its flexible leg supports, it allows individual fitting and can be positioned to accommodate the user's specific needs.

FEATURES

- **Tailored Support and Safety:** The RgoSling Toilet HighBack offers exceptional trunk support. It's designed with double belt buckles, ensuring secure attachment for additional safety.
- **Flexible Leg Supports:** The leg supports on the RgoSling Toilet HighBack are flexible, allowing for an individualized fitting and can be positioned in different ways to meet the user's needs.
- **Adaptable Sitting Angle:** The sling features several alternative loops on both the shoulder and leg straps. This design enables the adjustment of the sitting angle according to the user's needs, including allowing for an upright sitting position.
- **Comfort and Stability:** With its unique design, the RgoSling Toilet HighBack provides exceptional support in the sitting position. The adjustable belt buckles allow for optimal positioning of weight and pressure, ensuring maximum comfort, stability, and safety.

SPECIFICATIONS

Item Variant	Item number	Height (mm)	Recommended user weight (kg)	Waist size (mm)
XXS	1720600	600	12-17	400-650
XS	1720610	660	17-25	600-800

S	1720620	760	25-50	700-950
M	1720630	820	45-95	850-1100
L	1720640	890	90-165	1000-1250
XL	1720650	970	160-240	1200-1400
XXL	1720660	1030	230-300	1350-1600

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
 1300 633 468

support@medimarnt.com.au