



MOLIFT RGOSLING TOILET HIGHBACK PADDED

\$730.00 - \$1,040.00 Price range: \$730.00 through \$1,040.00



SKU: 17206-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Introducing the Molift RgoSling Toilet HighBack, a highback flexible sling engineered to provide optimal functions for toilet situations. This product is designed with large leg openings, making it easy to apply, remove, and facilitate toileting. With its flexible leg supports, it allows individual fitting and can be positioned to accommodate the user's specific needs.

FEATURES

- **Tailored Support and Safety:** The RgoSling Toilet HighBack offers exceptional trunk support. It's designed with double belt buckles, ensuring secure attachment for additional safety.
- **Flexible Leg Supports:** The leg supports on the RgoSling Toilet HighBack are flexible, allowing for an individualized fitting and can be positioned in different ways to meet the user's needs.
- **Adaptable Sitting Angle:** The sling features several alternative loops on both the shoulder and leg straps. This design enables the adjustment of the sitting angle according to the user's needs, including allowing for an upright sitting position.
- **Comfort and Stability:** With its unique design, the RgoSling Toilet HighBack provides exceptional support in the sitting position. The adjustable belt buckles allow for optimal positioning of weight and pressure, ensuring maximum comfort, stability, and safety.

SPECIFICATIONS

Item Variant	Item number	Height (mm)	Recommended user weight (kg)	Waist size (mm)
XXS	1720600	600	12-17	400-650
XS	1720610	660	17-25	600-800

S	1720620	760	25-50	700-950
M	1720630	820	45-95	850-1100
L	1720640	890	90-165	1000-1250
XL	1720650	970	160-240	1200-1400
XXL	1720660	1030	230-300	1350-1600