




## MOLIFT RGOSLING TOILET LOWBACK PADDED



**SKU:** 17205-H

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Introducing the Molift RgoSling Toilet LowBack, a flexible sling designed to provide optimal functionality for toilet situations. It features a large leg opening for toileting, making it easy to apply and remove while providing a user-friendly and individual-fitting experience.

### FEATURES

- **Trunk Support and Safety:** The RgoSling Toilet LowBack offers robust trunk support, reinforced by double belt buckles for extra security.
- **Adaptable Leg Supports:** The sling comes with flexible leg supports that facilitate individual fitting. These supports can be positioned in various ways to cater to the user's specific needs and function.
- **Adjustable Sitting Angle:** The RgoSling Toilet LowBack has numerous alternative loops on both the shoulder and leg straps. This feature allows adjustment of the sitting angle to meet the user's needs, including an upright sitting position.
- **Comfort and Stability:** With its unique design, the RgoSling Toilet LowBack provides outstanding support in the sitting position. The adjustable belt buckles allow for the optimal positioning of weight and pressure, maintaining excellent stability and safety for the user.

### SPECIFICATIONS

Item Variant	Item number	Height (mm)	Recommended user weight (kg)	Waist size (mm)
XXS	1720500	240	12-17	400-650
XS	1720510	260	17-25	600-800

S	1720520	280	25-50	700-950
M	1720530	300	45-95	850-1100
L	1720540	310	90-165	1000-1250
XL	1720550	320	160-240	1200-1400
XXL	1720560	330	230-300	1350-1600

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
 Brighton SA 5048  
 08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
 Ave, Elizabeth Vale SA 5112  
 08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
 Reynella SA 5161  
 08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
 03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
 03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
 1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)