



MOLIFT SMART 150 PORTABLE PATIENT LIFTER

\$7,400.00



SKU: M10200

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Awarded for Design Excellence in Europe, the Molift Smart 150 is a lightweight, portable patient lifter designed to be easily wheeled and transported. It is ideally suited for use in home care, for travelling, transporting in a car or to be easily stored.

With a lifting range of 27 - 168 cm (10.6" - 66.1") the Smart 150 provides an excellent maximum lifting height and at the same time gets low enough to facilitate lifting from the floor.

The 4-point sling suspension is standard and guarantees the patient hoist gives a comfortable lifting position for the user with ample space around the head area. This eliminates any risk of the head bumping into the suspension or the user swinging side to side.

With the combination of sling and suspension, the user will slide into a natural, correct seating position, from lying to sitting and vice versa.

The Smart 150 has a unique design with an excellent centre of balance and together with the small footprint makes it easy to manoeuvre even in tight spaces.

FEATURES

- Comfortable and safe transfer
- 4-point sling bar for more spacious lifting position
- Lifts easily from floor
- Small footprint makes it easy to manoeuvre in narrow spaces
- Low weight

- Folds and unfolds easily without the need for tools
- Easy and flexible to transport when travelling
- Built in software calculates when the lifter requires servicing
- Electric | manual lowering functions
- Emergency stop
- Includes: hand control | battery | charger | sling bar
- [Molift Hard Travel Case For Smart Lifter 150 available](#)

SPECIFICATIONS

Molift Smart 150

Item number	M10200
Battery	NiMH 14.4 V-2.6 Ah
Battery charge time (hours)	3
Hoisting height max (mm)	1665
Hoisting height min (mm)	210
Hoisting interval (mm)	1410
Hoisting speed 75 kg (mm/sec)	60
Leg height (mm)	110
Length (mm)	1200
Material	Aluminium, steel
Maximum user weight (kg)	150
Rated performance	40 hoists at: 75 kg and 500 mm
Service software	Molift Service Tool
Weight (kg)	26
Width chassis, legs open inner (mm)	925
Width chassis, outer (mm)	660