



MOLIFT SMART 150 PORTABLE PATIENT LIFTER



SKU: M10200

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Awarded for Design Excellence in Europe, the Molift Smart 150 is a lightweight, portable patient lifter designed to be easily wheeled and transported. It is ideally suited for use in home care, for travelling, transporting in a car or to be easily stored.

With a lifting range of 27 – 168 cm (10.6" – 66.1") the Smart 150 provides an excellent maximum lifting height and at the same time gets low enough to facilitate lifting from the floor.

The 4-point sling suspension is standard and guarantees the patient hoist gives a comfortable lifting position for the user with ample space around the head area. This eliminates any risk of the head bumping into the suspension or the user swinging side to side.

With the combination of sling and suspension, the user will slide into a natural, correct seating position, from lying to sitting and vice versa.

The Smart 150 has a unique design with an excellent centre of balance and together with the small footprint makes it easy to manoeuvre even in tight spaces.

FEATURES

- Comfortable and safe transfer
- 4-point sling bar for more spacious lifting position
- Lifts easily from floor
- Small footprint makes it easy to manoeuvre in narrow spaces
- Low weight

- Folds and unfolds easily without the need for tools
- Easy and flexible to transport when travelling
- Built in software calculates when the lifter requires servicing
- Electric | manual lowering functions
- Emergency stop
- Includes: hand control | battery | charger | sling bar
- [Molift Hard Travel Case For Smart Lifter 150 available](#)

SPECIFICATIONS

Molift Smart 150

Item number	M10200
Battery	NiMH 14.4 V-2.6 Ah
Battery charge time (hours)	3
Hoisting height max (mm)	1665
Hoisting height min (mm)	210
Hoisting interval (mm)	1410
Hoisting speed 75 kg (mm/sec)	60
Leg height (mm)	110
Length (mm)	1200
Material	Aluminium, steel
Maximum user weight (kg)	150
Rated performance	40 hoists at: 75 kg and 500 mm
Service software	Molift Service Tool
Weight (kg)	26
Width chassis, legs open inner (mm)	925
Width chassis, outer (mm)	660

For further information go to **[medimart.com.au](https://www.medimart.com.au)** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au