




MONTAR WALL MOUNTED CHANGE TABLE



SKU: N/A

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

Price on application. Please get in contact with us if you require more information or are wanting a quote.

MONTAR Wall-Mount Exam Couch The MONTAR fold-up exam couch will introduce space saving efficiencies and integrate seamlessly into your facilities.

The MONTAR incorporates Linak's state of the art electrics to achieve a 350mm low height and 950mm high height for ergonomic transfer of the patient from couch to stretcher. With industry grade upholstery and the frame is treated with an anti-bacterial polyester powder coating, the Montar Wall Mounted changetable is suitable for many applications such as

FEATURES

- First Aid
- Disabled Assistance
- Sports Injuries
- Medical Assessments
- Medical Treatment
- Examinations
- Training

The Montar is manufactured to exceed Australian Standards.